

This micro needling is a technique that builds collagen underneath the scars, not only this technique can be used for acne scarring but also treats fine lines and wrinkles, skin texture, pore size, brown spots, stretch marks, and pigment issues.

Nanoroller 196-needles gold plated: **Sizes available : 0.25 mm, 0.5 mm, 1.0 mm, 1.5 mm, 2.0 mm, 2.5 mm**

**NanoRoller** is a device that is used for collagen induction therapy, drug delivery, rejuvenates & repair scars and stretch marks.

#### Direction

- Wash & Degrease the target area
- Apply nano meso solution
- Sterile the nanoroller before use
- Glide the nanoroller gently on target area using the violin method
- Let nano meso solution penetrate, wait for 2 minutes followed by gentle massage

#### Warning

- Single-use only
- Keep away from children
- Stop using if irritation occurs
- Do not use it on irritated or infected skin
- Store in cool & dry place

0.25 Mm	● ○ ○ ○ ○
0.5 Mm	● ● ○ ○ ○
1.5 Mm	● ● ● ○ ○
2.5 Mm	● ● ● ● ●



**Hair Regrowth**

Max 0.5mm

**Deep Forehead wrinkle care**

Max 0.5mm

**Below the eyes**

Max 0.5mm  
dark circles, freckles

**The Bridge of the nose**

Max 0.5mm, large pores

**Philtrum**

Max 0.5mm, Elasticity

**Middle wrinkle of the forehead**

Max 1.5mm

**Temple**

Max 0.5mm

**Cheek**

Max 1.0mm, acne spots  
large pores, elasticity

**Fate Wrinkle (Deep Wrinkle)**

Max 1.0mm, elasticity

**Chin**

Max 0.5mm

LEADERMA  
Leader in Dermatology